|                             |    |
|-----------------------------|----|
| Indications for Use         | 2  |
| Indikationer for anvendelse | 3  |
| Gebruiksindicaties          | 4  |
| Indications d'emploi        | 5  |
| Anwendungsgebiete           | 6  |
| Ενδείξεις χρήσης            | 7  |
| Felhasználási javallatok    | 8  |
| Indicazioni per l'uso       | 9  |
| Bruksindikasjoner           | 10 |
| Indicações de utilização    | 11 |
| Indicaciones de uso         | 12 |
| Indikationer för användning | 13 |
| Kullanım Endikasyonları     | 14 |

---

**Posi-Stick® Introducer System**

---

**Indications for Use:**

The Posi-Stick® Introducer System is indicated for use in percutaneous introduction and placement of catheters and guidewires.

**Intended Patient Population and Group:**

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

**Device Description:**

The Posi-Stick® Introducer consists of an inner stiffening cannula, dilator, and an outer dilator. The stiffening cannula accepts a .018" guidewire, and is design to be inserted into the inner dilator. The device assembly may have a silicon or hydrophilic coating.

**Kit configurations may include the following devices:**

**Guidewire:** There are two configurations; Mandrel or Fully Coiled. Mandrel guidewire construction uses a tapered core wire which is soldered or welded to a coil at the distal end. Fully Coiled guidewire construction uses a tightly wound coil which surrounds a tapered core and a safety ribbon.

**Needle:** Introducer needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the target anatomy.

**Contraindications:**

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel or target anatomy. There is increased risk of pneumothorax for the patient who has severe chronic lung disease.

**Potential Complications:**

The potential complications related to the use of the introducer include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis hematoma formation, hemothorax, vessel erosion, trauma to vessels, inadvertent arterial or venous puncture, local or systemic infection/sepsis.

**Precautions:**

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

**USA CAUTION**

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

**Cautions:**

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Do not expose to organic solvents, eg. alcohol. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

**Warnings:**

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.
- Do not soak in hydrophilic coated Posi-Stick® Introducer System in alcohol.
- To ensure lubricity, do not leave the hydrophilic coated Posi-Stick® Introducer System hydrated for long periods of time.
- Avoid abrasion of the coating of the hydrophilic coated Posi-Stick® Introducer System.
- During insertion avoid contact with bone, cartilage and scar tissue, which can damage the needle, guidewire or introducer

**USE STERILE TECHNIQUE.** A suggested procedure:

1. Peel open package and place contents on sterile field. Inspect introducer and accessories for defects. Do not use any defective devices.
2. To remove air, flush the introducer with normal saline.  
**For Hydrophilic Coated Posi-Stick® Introducer System only:** soak the introducer with saline for 15 to 20 sec. in order to completely wet and hydrate the surface. When hydrated (wetted with saline or blood), the introducer is very lubricious. Use sterile gauze soaked with saline to help with keeping the guidewire hydrated and to aid in its handling.
3. Prep skin and drape in area of anticipated puncture site as desired.
4. Insert 21 ga trocar needle using sterile technique. The needle position should be verified.
5. The angle of the needle should be adjusted depending on the patient's build.
6. Remove trocar leaving needle cannula in place.
7. Advance the .018" guidewire through the introducer needle. Advance guidewire to required depth using ultrasound guidance, CT or fluoroscopy. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire location is suggested.
8. Hold guidewire in place and remove needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.  
**CAUTION:** Do not allow Guidewire to advance totally into patient.
9. Advance the Posi-Stick® Introducer over the .018" guidewire into the desired position.
10. Release and remove the stiffening cannula and the inner dilator, leaving the outer sheath and .018" guidewire in position.
11. Advance the .038" or .035" guidewire through the outer sheath (alongside the .018" guidewire) into the desired position.
12. Remove the .018" guidewire.
13. Remove the outer sheath leaving the .038" or .035" guidewire in place.
14. Utilize guidewire for placement of further diagnostic or interventional devices.

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
Posi-Stick<sup>®</sup> Introducer System

|  |   |  |   |  |
|--|---|--|---|--|
| <p><b>Rx ONLY</b></p> <p>Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)</p> |  <p>Do Not Use if Package Is Damaged</p> |  <p>Single Use Only</p> |  <p>Not Made With Natural Rubber Latex</p> |  <p>Do Not Re-Sterilize</p>   |
|  <p>Ethylene Oxide Sterilized</p>                         |  <p>Instructions for Use</p>             |  <p>Caution</p>         |  <p>Non-Pyrogenic</p>                      |  <p>Store In A Dark Place</p>   |
|  <p>Keep Dry</p>  |  <p>Part Number</p>                      |  <p>Lot Number</p>      |  <p>Manufacturing Date</p>                 |  <p>Use By Date</p>   |
|  <p>Manufacturer</p>                                      |  <p>Authorized Representative</p>        |  |  <p>Quantity</p>                           |  <p>Medical Device</p>  |
| <p>NiTi</p> <p>Nitinol</p>   | <p>NiTi [Au]</p> <p>Nitinol Gold</p>  | <p>NiTi [Pt]</p> <p>Nitinol Platinum</p>   |  <p>Single Sterile Barrier system</p>      |  <p>Single Sterile Barrier System with Protective Packaging Outside</p> |