		
EN Catalog number CS Katalogové číslo DA Katalognummer DE Katalognummer EL Αριθμός καταλόγου ES Número de catálogo FR Numéro de catalogue IT Numero di catalogo NL Catalogusnummer NO Katalognummer PT Referência SK Katalógové číslo SV Katalognummer TR Katalog numarası ZH 产品目录号	EN Date of manufacture CS Datum výroby DA Fremstillingsdato DE Herstellungsdatum EL Ημερομηνία κατασκευής ES Fecha de fabricación FR Date de fabrication IT Data di produzione NL Productiedatum NO Produksjonsdato PT Data de fabrico SK Datum výroby TR Tilverkningsdatum SV Utlätim tärhi ZH 生产日期	EN Batch code CS Kód šarže DA Batchkode DE Chargencode EL Κωδικός παρτίδας ES Código de lote FR Code du lot IT Codice lotto NL Batchcode NO Partikode PT Código do lote SK Kód šarže SV Batchkod TR Parti kodu ZH 批号

		
EN Telephone Number CS Telefonní číslo DA Telefonnummer DE Telefonnummer EL Αριθμός τηλεφώνου ES Número de teléfono FR Numéro de téléphone IT Numero di telefono NL Telefoonnummer NO Telefonnummer PT Número de telefone SK Telefónne číslo SV Telefonnummer TR Telefon Numarası ZH 电话号码	EN Fax Number CS Číslo faxu DA Faxnummer DE Faxnummer EL Αριθμός φαξ ES Número de fax FR Numéro de fax IT Numero di fax NL Faxnummer NO Faksnummer PT Número de fax SK Faxové číslo SV Faxnummer TR Faks Numarası ZH 传真号码	EN Use-by date CS CS Datum použitelnosti DA Sidste anvendelsesdato DE Verfallsdatum EL Ημερομηνία λήξης ES Fecha de caducidad FR Date de péremption IT Data di scadenza NL Houdbaarheidsdatum NO Bruk før-dato PT Prazo de validade SK Dátum spotreby SV Sista förbrukningsdag TR Son kullanma tarihi ZH 有效期

		UreSil, LLC, 5418 W. Touhy Ave, Skokie, IL 60077 U.S.A.
EN Unique Device Identifier CS Jediný identifikační prostředek DA Unik enhedsidentifikator DE Einzigartige Geräteerkennung EL Μοναδικό αναγνωριστικό συσκευής ES Número único de identificación de producto FR Identifiant unique du dispositif IT Identificativo unico del dispositivo NL Uniek apparaat-ID NO Unik enhetsidentifikator PT Identificação Única do Dispositivo SK Unikátny identifikátor pomôcky SV Unik identifikationskod TR Benzersiz Cihaz Tanımlayıcısı ZH 设备唯一标识符	EN Procedure Pack Producer CS Postup výrobce obalu DA Producent af procedurepakke DE Hersteller von Behandlungseinheiten EL Παραγωγός Πακέτου Διαδικασίας ES Productor del conjunto de procedimiento FR Producteur de packs de procédures IT Produttore del pacco procedurale NL Producent procedurepakketten NO Prosedyrepakkens produsent PT Procedimento Embalagem Produtor SK Výrobca súpravy pomôck SV Tillverkare procedurförpackning TR Prosedür Paketleri Üreticisi ZH 工序包装生产商	

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Fax

Distribútor

**Posi-Stick® Introducer System****Indications for Use:**

The Posi-Stick® Introducer System is indicated for use in percutaneous introduction and placement of catheters and guidewires.

Intended Patient Population and Group:

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

Device Description:

The Posi-Stick® Introducer consists of an inner stiffening cannula, dilator, and an outer dilator. The stiffening cannula accepts a .018" guidewire, and is design to be inserted into the inner dilator. The device assembly may have a silicon or hydrophilic coating.

Kit configurations may include the following devices:

Guidewire: There are two configurations; Mandrel or Fully Coiled. Mandrel guidewire construction uses a tapered core wire which is soldered or welded to a coil at the distal end. Fully Coiled guidewire construction uses a tightly wound coil which surrounds a tapered core and a safety ribbon.

Needle: Introducer needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the target anatomy.

Contraindications:

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel or target anatomy. There is increased risk of pneumothorax for the patient who has severe chronic lung disease.

Potential Complications:

The potential complications related to the use of the introducer include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis hematoma formation, hemothorax, vessel erosion, trauma to vessels, inadvertent arterial or venous puncture, local or systemic infection/sepsis.

Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

Cautions:

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, deter mine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Do not expose to organic solvents, eg. alcohol. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

Warnings:

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.
- Do not soak in hydrophilic coated Posi-Stick® Introducer System in alcohol.
- To ensure lubricity, do not leave the hydrophilic coated Posi-Stick® Introducer System hydrated for long periods of time.
- Avoid abrasion of the coating of the hydrophilic coated Posi-Stick® Introducer System.
- During insertion avoid contact with bone, cartilage and scar tissue, which can damage the needle, guidewire or introducer

USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field. Inspect introducer and accessories for defects. Do not use any defective devices.
2. To remove air, flush the introducer with normal saline.
For Hydrophilic Coated Posi-Stick® Introducer System only: soak the introducer with saline for 15 to 20 sec. in order to completely wet and hydrate the surface. When hydrated (wetted with saline or blood), the introducer is very lubricious. Use sterile gauze soaked with saline to help with keeping the guidewire hydrated and to aid in its handling.
3. Prep skin and drape in area of anticipated puncture site as desired.
4. Insert 21ga trocar needle using sterile technique. The needle position should be verified.
5. The angle of the needle should be adjusted depending on the patient's build.
6. Remove trocar leaving needle cannula in place.
7. Advance the .018" guidewire through the introducer needle. Advance guidewire to required depth using ultrasound guidance, CT or fluoroscopy. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire location is suggested.
8. Hold guidewire in place and remove needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
CAUTION: Do not allow Guidewire to advance totally into patient.
9. Advance the Posi-Stick® Introducer over the .018" guidewire into the desired position.
10. Release and remove the stiffening cannula and the inner dilator, leaving the outer sheath and .018" guidewire in position.
11. Advance the .038" or .035" guidewire through the outer sheath (alongside the .018" guidewire) into the desired position.
12. Remove the .018" guidewire.
13. Remove the outer sheath leaving the .038" or .035" guidewire in place.
14. Utilize guidewire for placement of further diagnostic or interventional devices.

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
 Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)	 Do Not Use If Package Is Damaged	 Single Use Only	 Not Made With Natural Rubber Latex	 Do Not Re-Sterilize
 Ethylene Oxide Sterilized	 Instructions for Use	 Caution	 Non-Pyrogenic	 Store In A Dark Place
 Keep Dry	 Part Number	 Lot Number	 Manufacturing Date	 Use By Date
 Manufacturer	 Authorized Representative	 Quantity	 Medical Device	 MR Unsafe
 Nitinol	 NiTi [Au] Nitinol Gold	 NiTi [Pt] Nitinol Platinum	 Single Sterile Barrier system	 Single Sterile Barrier System with Protective Packaging Outside

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 Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)	 Do Not Use If Package Is Damaged	 Single Use Only	 Not Made With Natural Rubber Latex	 Do Not Re-Sterilize
 Ethylene Oxide Sterilized	 Instructions for Use	 Caution	 Non-Pyrogenic	 Store In A Dark Place
 Keep Dry	 Part Number	 Lot Number	 Manufacturing Date	 Use By Date
 Manufacturer	 Authorized Representative	 Quantity	 Medical Device	 MR Unsafe
NiTi Nitinol	NiTi [Au] Nitinol Gold	NiTi [Pt] Nitinol Platinum	Single Sterile Barrier system	Single Sterile Barrier System with Protective Packaging Outside

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DILATOR

Indications for Use:

These dilators are used for the percutaneous introduction of guidewires into the peripheral vasculature.

Device Description:

The dilator introducer consists of a tipped radiopaque tube with a molded hub that accepts up to a .038 inch. guidewires

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Contraindications:

Use of the dilator is contraindicated if the patient has a known or suspected obstruction in the vessel. There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

Potential Complications:

The potential complications related to the use of the dilator include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis, hematoma formation, hemothorax, vessel erosion, trauma to vessels, sepsis..

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Cautions:

- This procedure should only be performed by professionals thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire through the needle, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

USE STERILE TECHNIQUE.

A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture site as desired.
3. Distend the vessel of choice following standard hospital practice for venipuncture. The vein will be much easier to locate if the patient is well hydrated.
4. Insert needle into vessel. The needle position should be verified by observing venous blood return.
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
6. Aspirate the puncture needle using the syringe.
7. Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire position is suggested.
8. Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.
Caution: If the guidewire must be withdrawn while the needle is inserted, remove both the needle and wire as a unit to prevent the needle from damaging or shearing the guidewire.
9. Thread the dilator over the guidewire.
10. Advance the dilator with a twisting motion over the guidewire and into the vessel, to dilate tissue in vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
11. Slowly withdraw and remove the dilator, while holding the guidewire in position and insert catheter or over guidewire.
12. Remove guidewire leaving catheter in position.

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