DA

Indføringskanyle til UreSil THORA-VENT Thoraxventili 11F, 13 cm

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
MD

Insertion Cannula for UreSil®
THORA-VENT® Thoracic Vent 11F x 13cm

INTENDED USE: This over-the-wire cannula is intended for Seldinger Introduction of UreSil THORA-VENT - 11F x 13cm. This product is intended for transient use (less than 24 hours).

INTENDED USERS: This product is intended only for use by trained clinicians.

Insertion Cannula

INSTRUCTIONS FOR USE:

This package contains instructions for the specific use of the over-the-wire cannula. For additional instructions related to the use of the Thoracic Vent, see the instructions for use that are included in the Thoracic Vent Procedure Tray.

- Select and prepare the site as described in the instructions provided in the Thoracic Vent Procedure Tray.
- Using the scalpel, make a small incision at the selected site.
- Using Seldinger technique, dilate the tract to approximately one French size larger than the catheter on the Thoracic Vent.
- Remove the protective tube from the cannula and insert it into the self-sealing port. Lock the cannula into the Thoracic Vent.
- Peel away the center portions of the paper cover from the adhesive patch on the Thoracic Vent.
- Advance the distal end of the catheter/cannula assembly into the Thoracic Vent over the guidewire and into the pleural space. When the tip of the catheter enters the pleural space, stop advancing, unlock the cannula and introduce the full length of the catheter. Remove the cannula and the guidewire. Do not re-insert the cannula until the Thoracic Vent has been completely removed from the patient.
- Read the instructions for use in the Thoracic Vent Procedure Tray for additional information.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil and THORA-VENT are registered trademarks of UreSil, LLC, U.S. Patent no. 6,770,063.

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