DA

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
MD

UreSil ® "Y" Connector Tubing Set

INTENDED USE: For connecting a percutaneous drainage catheter to a drainage bag or system. For bilateral drainage with a gravity drainage or suction drainage system. This product is intended for short-term use (less than 30 days).

INTENDED PATIENT POPULATION: This product is intended to be used with patients with placed drainage catheters.

INTENDED USERS: This product is intended only for use by trained clinicians.


"Y" Connector Tubing Set

CLINICAL BENEFITS: Work as an accessory in conjunction with a percutaneous drainage catheter and a Gravity Drainage bag or system to remove fluids using gravity or in conjunction with a percutaneous drainage catheter and a Suction Drainage bag or system to remove fluids using active suction or gravity.

PERFORMANCE CHARACTERISTIC: The TCY2 "Y" Connector Tubing Set is composed of an inlet tube with a luer connector and a "Y" connector and is used to allow bilateral percutaneous catheter drainage into a single Gravity Drainage or Suction Drainage bag or system.

INTENDED FOR USE WITH: This accessory for a medical device is intended for use with bilateral percutaneous drainage catheters with standard luer lock and a Gravity Drainage or Suction Drainage bag or system.

INSTRUCTIONS FOR USE:

CONTRAINDICATIONS: Currently no specific contraindications are known for this accessory for a medical device.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction. Do not use if the sterile package is damaged or is unintentionally opened before use.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil is a registered trademark of UreSil, LLC.

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