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

INTENDED USE: For collecting drainage fluids. The Gravity Drainage Bags are intended to be used to collect fluids drained from abscesses, the hepatic system, and the kidneys. This product is intended for short-term use (less than 30 days). Replace bag every 5 to 7 days or sooner if clinically indicated, for example, if malodorous or damaged.

INTENDED PATIENT POPULATION: This product is intended to be used with patients with placed drainage catheters.

INTENDED USERS: This product is intended only for use as directed by a trained clinician. A nurse, caregiver, or patient may drain the fluids using the provided empty port as directed by a trained clinician.

NOTE: Graduation marks on bag not intended to measure – for reference only.

CLINICAL BENEFITS: The Gravity Drainage Bags are available in a variety of sizes and configurations. They can attach to any standard luer locking drainage catheter and can remove fluids using gravity. Additional clinical benefits:

- Simple and intuitive to operate.
• Compact and self-contained system allows full patient ambulation.
• During activation, body fluids are not aerosolized because all fluids are contained in the system.
• Closed system design reduces the chances of cross-contamination.
• Anti-reflux valve prevents liquids and air from backing up into the catheter.

PERFORMANCE CHARACTERISTIC: Gravity Drainage Bags have an inlet tube with a luer lock fitting for connection to a standard luer lock drainage catheter, anti-reflux valve, empty port and Velcro belt system. Flow to the bag occurs via gravity drainage as long as the bag is below the level of the drainage site.

ACCESSORIES TO BE USED WITH THE DEVICE: The PDC-36 Percutaneous Drainage Catheter Connector Tubing is comprised of an inlet tube with a luer connector and is used for connecting a percutaneous catheter to a drainage bag or system. The TCY2 "Y" Connector Tubing Set is comprised of an inlet tube with a luer connector and a "Y" connector and is used for bilateral drainage with a gravity-drainage or suction-drainage system. UFA-LL is comprised of a tapered funnel and a female luer connector and is used to connect UreSil Uti-Ilo Large Bore non-locking catheters (REF LBNL-series) to a standard luer.

INSTRUCTIONS FOR USE:

1. Close the port prior to use.
2. Thread the belt as shown.
3. Connect the drainage catheter to the bag tubing as shown.
4. If altering the bag tube length, using sterile technique, remove the luer lock, cut the tubing and re-attach the luer lock.
5. Suspend the bag at a location below the drainage site as shown.
6. Replace bag every 5 to 7 days or sooner if clinically indicated, for example, if malodorous or damaged.

- CAUTION:
1. Use care if emptying bag to avoid contact with drained fluids.
2. Consult a healthcare professional if no fluid is draining into the bag or if fluid is leaking from the bag.

CONTRAINDICATIONS: Currently no specific contraindications are known for the Gravity Drainage Bags.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction. Sterile if package is unopened and undamaged. Do not use if the sterile package is damaged or is unintentionally opened before use.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

Optimus and UreSil are registered trademarks of UreSil, LLC.

8-19-200 rev. 04/10/2023


Table with 3 columns: REF, LOT, and CAUTION symbols. Each column contains multilingual text for reference, lot, and caution information.

Table with 2 columns: USA and EN. Each column contains multilingual text for country of manufacture and date of manufacture information.

Table with 2 columns: EN and CS. Each column contains multilingual text for quantity and date of manufacture information.

Table with 2 columns: EN and CS. Each column contains multilingual text for latex-free status information.

Table with 2 columns: EC REP and UDI. Each column contains multilingual text for representative and unique device identifier information.

Table with 2 columns: EN and CS. Each column contains multilingual text for importer information.

Table with 2 columns: EN and CS. Each column contains multilingual text for sterilization instructions.

Table with 2 columns: EN and CS. Each column contains multilingual text for package damage and use instructions.

Table with 2 columns: EN and CS. Each column contains multilingual text for latex-free status information.

Table with 2 columns: EN and CS. Each column contains multilingual text for manufacturer and UDI information.

Table with 2 columns: EN and CS. Each column contains multilingual text for telephone and fax contact information.

Table with 2 columns: EN and CS. Each column contains multilingual text for sterilization and reuse instructions.

Table with 2 columns: EN and CS. Each column contains multilingual text for sterilization and reuse instructions.

Table with 2 columns: EN and CS. Each column contains multilingual text for working length information.

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