
See the pouch label on this product for specific contents.

INTENDED USE: This product is intended for short-term use (less than 30 days).

- General-purpose (GP/Mini-pigtail)/NEPHRO-CATH Nephrostomy: Percutaneous drainage of abscess, nephrostomy and other fluids.
- CHOLE-CATH Biliary: Percutaneous transhepatic biliary drainage.


INTENDED USERS: This product is intended only for use by trained clinicians.

LOCKING CATHETER: GP, NEPHROSTOMY AND BILIARY CATHETERS WITH LOCKING PIGTAIL

HARDWARE:

METAL STIFFENING CANNULA


SEPTUM CAP (IF PRESENT)

FLEXIBLE CANNULA (IF PRESENT)

TROCAR (IF PRESENT)

OBTURATOR (IF PRESENT)

NOTE: Depth marks (if present) not intended to measure – for reference only

INSTRUCTIONS FOR USE:

1. To maximize advantages of the surface, WET THE CATHETER before use.
2. Straighten the catheter when inserting the cannula (and trocar if appropriate). For locking catheters, pull the suture to remove slack before inserting the cannula.
3. Insert the catheter using the over-the-guidewire technique. Tract over-dilation will ease placement. Catheters with a trocar included can also be inserted using direct puncture.
4. For locking catheters, rotate the catheter clockwise while removing the cannula to form the pigtail.
5. When the catheter is properly positioned, remove the hardware.
6. For locking pigtail catheters, pull the suture to close the pigtail. Depending upon whether or not suture storage is desired, lock the pigtail as follows:
 - Option A: Locking the pigtail without storing suture: Move the locking slide proximally until an audible “click” is heard. The suture is now locked. Trim away excess suture at the exit point at the proximal end of the hub. Do not pull suture out of catheter until removing catheter from patient, as this may cause a leak.
 - Option B: Locking the pigtail with suture storage: Wind approximately 15 cm of suture into the suture channel for later use when removing or repositioning the catheter. To wrap the suture, pull it straight and then hook it in the notch at the bottom of the suture channel. After making the last wind, wrap the suture around the suture retention post. Move the locking slide proximally until an audible “click” is heard. The suture is now locked.
7. To remove the catheter, unlock the suture retention mechanism. To unlock the mechanism, depress the tab using an instrument and move the slide distally. If no suture was stored, the suture will now be released. If suture was stored, unwind the suture to complete the release. Cut the suture at the exit point at the proximal end of the hub and pull the catheter out of the patient. In order to maintain control of the suture, do not cut the catheter to release the pigtail.
8. For biliary and nephrostomy catheters, internal drainage can be achieved by placing the septum cap on the catheter hub. The septum cap contains a self-sealing needle puncture port for flushing.

CAUTION:

1. When engaging the catheter luer lock, do not overtighten. Do not use instruments to tighten connections. Overtightening the luer can result in luer crack which may require catheter replacement.
2. Dispose of sharps using approved sharps container in accordance with applicable regulations and institutional policy.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

URESIL, MINI-PIGTAIL, CHOLE-CATH, and Origin are registered trademarks of UreSil, LLC. GP and NEPHRO-CATH are trademarks of UreSil, LLC. U.S. patent no's. 8,177,773; 9,079,006B1.

REF

EN Catalog number
 CS Katalogové číslo
 DA Katalognummer
 DE Katalognummer
 EL Αριθμός καταλόγου
 ES Número de catálogo
 FI Luettelonumero
 FR Numéro de catalogue
 IT Numero di catalogo
 NL Catalogusnummer
 NO Katalognummer
 PT Referência
 SK Katalógové číslo
 SL Kataloška številka
 SV Katalognummer
 TR Katalog numarası
 ZH 产品目录号

LOT

EN Batch code
 CS Kód šarže
 DA Batchkode
 DE Chargencode
 EL Κωδικός παρτίδας
 ES Código de lote
 FI Eräkkoodi
 FR Code du lot
 IT Codice lotto
 NL Batchcode
 NO Partikode
 PT Código do lote
 SK Kód šarže
 SL Koda serije
 SV Batchkod
 TR Parti kodu
 ZH 批号

EN Caution
 CS Upozornění
 DA Forsigtig
 DE Achtung
 EL Προσοχή
 ES Precaución
 FI Huomio
 FR Avertissement
 IT Attenzione
 NL Opgelet
 NO Advarsel
 PT Atenção
 SK Upozornenie
 SL Pozor
 SV Varning
 TR Dikkat
 ZH 注意

EN Country of manufacture/ Date of manufacture
 CS Země výroby/ Datum výroby
 DA Fremstillingsland/ Fremstillingsdato
 DE Herstellungsland/ Herstellungsdatum
 EL Χώρα κατασκευής/ Ημερομηνία κατασκευής
 ES País de fabricación/ Fecha de fabricación
 FI Valmistusmaa/ Valmistuspäivä
 FR Pays de fabrication/ Date de fabrication
 IT Paese di produzione/ Data di produzione
 NL Land van productie/ Productiedatum
 NO Produksjonsland/ Produksjonsdato
 PT País de fabric/ Data de fabrico
 SK Krajina výroby/ Dátum výroby
 SL Država proizvodnje/ Datum proizvodnje
 TR Ambalaj hasarlıysa ürünü kullanmayın ve kullanım talimatlarına başvurun/ Tillverkningsdatum
 SV Tillverkningsland/ Üretim tarihi
 ZH 产地/ 生产日期

EN Quantity
 CS Množství
 DA Antal
 DE Menge
 EL Ποσότητα
 ES Cantidad
 FI Määrä
 FR Quantité
 IT Quantità
 NL Hoeveelheid
 NO Antall
 PT Quantidade
 SK Množstvo
 SL Količina
 SV Kvantitet
 TR Miktar
 ZH 数量

EN Not made with natural rubber latex
 CS Není vyrobeno z latexu z přírodního kaučuku
 DA Ikke fremstillet med naturgummilætex
 DE Nicht aus Naturkautschuklatex hergestellt
 EL Δεν κατασκευάζεται με λατέξ από φυσικό καουτσούκ
 ES No fabricado con látex de caucho natural
 FI Ei ole valmistettu luonnonkumilateksista
 FR Non fabriqué avec du latex de caoutchouc naturel
 IT Non realizzato con lattice di gomma naturale
 NL Niet gemaakt van natuurlijk rubberlatex
 NO Skal ikke brukes med naturlig gummilætex
 PT Não fabricado com látex de borracha natural
 SK Neobsahuje prírodný latex
 SL Ni narejeno iz naravnega kavčuka
 SV Inte tillverkad med naturlig gummilætex
 TR Doğal kauçuk lateksten yapılmamıştır
 ZH 并非天然橡胶材质

EN Do not re-use
 CS Nepoužívejte opakovaně
 DA Må ikke genbruges
 DE Nicht wiederverwenden
 EL Να μην επαναχρησιμοποιηθεί
 ES No reutilizar
 FI Älä käytä uudelleen
 FR Ne pas réutiliser
 IT Non riutilizzare
 NL Niet hergebruiken
 NO Ikke til gjenbruk
 PT Não reutilizar
 SK Nepoužívajte opakovane
 SL Ne uporabite znova
 SV Får inte återanvändas
 TR Yeniden kullanmayın
 ZH 不可重复使用

1-800-538-7374 (U.S.A.)

1-847-982-0106

EN Telephone Number
 CS Telefonní číslo
 DA Telefonnummer
 DE Telefonnummer
 EL Αριθμός τηλεφώνου
 ES Número de teléfono
 FI Puhelinnumero
 FR Numéro de téléphone
 IT Numero di telefono
 NL Telefoonnummer
 NO Telefonnummer
 PT Número de telefone
 SK Telefónne číslo
 SL Telefonska številka
 SV Telefonnummer
 TR Telefon Numarası
 ZH 电话号码

EN Fax Number
 CS Číslo faxu
 DA Faxnummer
 DE Faxnummer
 EL Αριθμός φαξ
 ES Número de fax
 FI Faksinumero
 FR Numéro de fax
 IT Numero di fax
 NL Faxnummer
 NO Faksnummer
 PT Número de fax
 SK Faxové číslo
 SL Številka faksa
 SV Faxnummer
 TR Faks Numarası
 ZH 传真号码

EN Do not re-sterilize
 CS Neresterilizujte
 DA Må ikke reesteriliseres
 DE Nicht erneut sterilisieren
 EL Να μην αποστειρωθεί εκ νέου
 ES No reesterilizar
 FI Älä steriloï uudelleen
 FR Ne pas restériliser
 IT Non risterilizzare
 NL Niet opnieuw steriliseren
 NO Skal ikke reesteriliseres
 PT Não reesterilizar
 SK Nesterilizujte
 SL Ne sterilizirajte znova
 SV Får inte omsteriliseras
 TR Yeniden sterilize etmeyin
 ZH 不可重复消毒

CH REP

MedEnvoy Switzerland, Gotthardstrasse 28,
 6302 Zug Switzerland

EN Authorised Representative in Switzerland/Importer in Switzerland
 CS Oprávněný zástupce ve Švýcarsku / Dovozece ve Švýcarsku
 DA Autoriseret repræsentant i Schweiz/importør i Schweiz
 DE Bevollmächtigter Vertreter in der Schweiz / Importeur in der Schweiz
 EL Εξουσιοδοτημένος εκπρόσωπος στην Ελλάδα/ εισαγωγέας στην Ελλάδα
 ES Representante autorizado en Suiza/Importador en Suiza
 FI Valtuutettu edustaja Sveitsissä / Maahantuojia Sveitsissä.
 FR Représentant autorisé en Suisse / Importateur en Suisse
 IT Rappresentante autorizzato in Svizzera/Importatore in Svizzera
 NL Geautoriseerde vertegenwoordiger in Zwitserland/importeur in Zwitserland
 NO Autorisert representant i Sveits / importør i Sveits
 PT Representante autorizado na Suíça/ Importador na Suíça
 SK Autorizovaný zástupca vo Švajčiarsku/dovožca vo Švajčiarsku
 SL Pooblaščen zastopnik v Švici/Uvoznik v Švici
 SV Auktoriserad representant i Schweiz/importör i Schweiz
 TR İsviçre'deki Yetkili Temsilci/ İsviçre'deki İthalatçı
 ZH 瑞士授权代表/瑞士进口商

EN Do not use if package is damaged and consult instructions for use
CS Nepoužívejte, pokud je obal poškozený, a přečtěte si návod k použití
DA Må ikke bruges, hvis pakken er beskadiget, og se brugsanvisningen
DE Nicht verwenden, wenn die Verpackung beschädigt ist, und Gebrauchsanweisung beachten
EL Μην χρησιμοποιείτε το προϊόν εάν η συσκευασία έχει υποστεί ζημιά και συμβουλευτείτε τις οδηγίες χρήσης
ES No utilizar si el envase está dañado y consultar las instrucciones de uso
FI Älä käytä, jos pakkaus on vahingoittunut, ja katso käyttöohjeet
FR Ne pas utiliser si l'emballage est endommagé et consulter le mode d'emploi
IT Non utilizzare se la confezione è danneggiata e consultare le istruzioni per l'uso
NL Niet gebruiken als de verpakking beschadigd is. Raadpleeg de gebruiksaanwijzing
NO Ikke bruk hvis pakke er skadet, og les bruksanvisningen
PT Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização
SK Nepoužívať, ak je obal poškodený, a prečítajte si návod na použitie
SL Ne uporabljajte, če je embalaža poškodovana, in glejte navodila za uporabo
SV Använd inte produkten om förpackningen är skadad och se bruksanvisningen
TR Ambalaj hasarlıysa ürünü kullanmayın ve kullanım talimatlarını başvurun
ZH 包装如有损坏不可使用，请参考使用说明书

EN Sterilized using ethylene oxide/ Single sterile barrier system
CS Sterilizováno ethylenoxidem/ Jednoduchý sterilní bariérový systém
DA Steriliseret ved anvendelse af ethylenoxid/ Enkelt steril barrieresystem
DE Mit Ethylenoxid sterilisiert/ Einfaches Sterilbarrieresystem
EL Αποστειρώθηκε με χρήση αιθυλενοξειδίου/ Σύστημα μονού αποστειρωμένου φραγμού
ES Esterilizado con óxido de etileno/ Sistema de barrera estéril simple
FI Steriloitu etyleenioksidilla/ Yksittäinen steriili estojärjestelmä
FR Stérilisé à l'oxyde d'éthylène/ Système de barrière stérile unique
IT Sterilizzato con ossido di etilene/ Sistema di barriera sterile singola
NL Gesteriliseerd met ethyleenoxide/ Enkelvoudig steriel barrieresysteem
NO Sterilisert med etylenoksid/ System med enkelt steril barriere
PT Esterilizado por óxido de etileno/ Sistema de barreira estéril simples
SK Sterilizované etylénoxidom/ Systém jednej sterilnej bariéry
SL Sterilizirano z etilen oksidom/ Enojni sterilni sistem pregrad
SV Steriliserad med etylenoxid/ Sterilt barriärsystem med ett lager
TR Etilen oksit kullanılarak sterilize edilmiştir/ Tek steril bariyerli sistem
ZH 使用环氧乙烷消毒/ 单一无菌屏障系统

EN Locking pigtail
CS Zajišťující pigtail
DA Låsende grisehale
DE Pigtail-Verriegelung
EL Άκρο τύπου pigtail με δυνατότητα κλειδώματος
ES Pigtail de bloqueo
FI Lukittava letkukatetri
FR Verrouillage en queue de cochon
IT Bloccaggio "a coda di maiale" (pigtail)
NL Pigtail met vergrendeling
NO Låsende grisehale
PT Pigtail de retenção
SK Uzamykací pigtail
SL Zaklepni jeziček
SV Låsande pigtail-koppling
TR Kilitli pigtail
ZH 锁定猪尾管

EN Working length
CS Pracovní délka
DA Arbejds længde
DE Arbeitslänge
EL Λειτουργικό μήκος
ES Longitud de trabajo
FI Työskentelypituus
FR Longueur de travail
IT Lunghezza di lavoro
NL Werklengte
NO Arbeidslengde
PT Comprimento de trabalho
SK Pracovní dĺžka
SL Delovna dolžina
SV Arbetslängd
TR Çalışma uzunluğu
ZH 工作长度

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