MD

UreSil


EN

UreSil® POSI-STICK® 2 One-Step Introducer Set

INTENDED USE: This device is intended for placement of a percutaneous drainage catheter. This product is intended for transient use (less than 24 hours).

INTENDED PATIENT POPULATION: This product is intended to be used with patients with placed drainage catheters.

INTENDED USERS: This product is intended only for use by trained clinicians.

CLINICAL BENEFITS: Insertion of a percutaneous drainage catheter using the over-the-guidewire technique.

PERFORMANCE CHARACTERISTIC: The Posi-Stick One-Step Introducer Set includes a trocar, dilator, and guidewire that allows the clinician to insert a percutaneous drainage catheter using the over-the-guidewire technique.

TO BE USED WITH: The Posi-Stick One-Step Introducer Set (PS2-1S) is intended to be used with percutaneous drainage catheters.

INSTRUCTIONS FOR USE:

1. Insert the trocar with the dilator into the drainage site using standard technique.
2. Remove the trocar leaving the flexible dilator in place.
3. Advance the guidewire through the dilator into the drainage site.
4. Remove the dilator leaving the guidewire in place.

CAUTION:

1. Withdrawal of guidewires should be smooth and without resistance. If resistance is felt, carefully remove wire(s) and system as a unit.
2. Avoid bending, kinking, and uses other than percutaneous catheter placement.
3. Guidewires should be routinely inspected prior to use. Should any deformities be present, discard. Do not use deformed guidewires.
4. Should resistance occur during insertion, DO NOT advance the guidewire.
5. Avoid withdrawing a guidewire through a metal needle. The teflon coating can shear against the needle bevel.

CONTRAINDICATIONS: Currently no specific contraindications are known for the Posi-Stick One-Step Introducer Set.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction. Sterile if package is unopened and undamaged. Do not use if the sterile package is damaged or is intentionally opened before use.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil and Posi-Stick are registered trademarks of UreSil, LLC.

EN Do not resterilize
ES No reesterilizar
FR Ne pas restériliser
PT Não reesterilizar
ZH 不可重复消毒

EN Do not re-use
ES No reutilizar
FR Ne pas réutiliser
PT Não reutilizar
ZH 不可重复使用

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