MD

UreSil®

EN

UreSil® Percutaneous Drainage Catheter
Connector Tubing 36" (90cm)

INTENDED USE: For connecting a percutaneous drainage catheter to a drainage bag or system. This product is intended for short-term use (less than 30 days).

INTENDED PATIENT POPULATION: This product is intended to be used with patients with placed drainage catheters.

INTENDED USERS: This product is intended only for use by trained clinicians.

Percutaneous Drainage Catheter Connector Tubing


CLINICAL BENEFITS: Work as an accessory in conjunction with a percutaneous drainage catheter and a Gravity Drainage bag or system to remove fluids using gravity or in conjunction with a percutaneous drainage catheter and a Suction Drainage bag or system to remove fluids using active suction or gravity.

PERFORMANCE CHARACTERISTIC: The PDC-36 Percutaneous Drainage Catheter Connector Tubing is comprised of an inlet tube with a luer connector and is used for connecting a percutaneous drainage catheter to a drainage bag or system with a tapered connector.

INTENDED FOR USE WITH: This accessory for a medical device is intended for use with a percutaneous drainage catheter with standard luer lock and a Gravity Drainage or Suction Drainage bag or system with a tapered connector.

INSTRUCTIONS FOR USE:

1. Cut the flared end to fit tapered connector or drainage bag or fluid removal system.
2. Attach the luer lock connector to the catheter hub.

CONTRAINDICATIONS: Currently no specific contraindications are known for this accessory for a medical device.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction. Sterile if package is unopened and undamaged. Do not use if the sterile package is damaged or is unintentionally opened before use.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil is a registered trademark of UreSil, LLC.

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EC REP

EMERGO EUROPE, Westervoortsedijk 60, 6827
AT Arnhem, The Netherlands

EN Authorized representative in the European Community / European Union

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EN Importer in the European Community / European Union

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