
## **UreSil Catheter Sales Kit**

Contains:

IFU for Posi-Stick Introducer  
System LAB-1032-00, Rev. C

IFU for Dilator Pack LAB-1117-00,  
Rev. B

MedEnvoy Switzerland, Gotthardstrasse 28,  
6302 Zug Switzerland

EN Authorised Representative in Switzerland/Importer in Switzerland  
CS Oprávněný zástupce ve Švýcarsku / Dovozece ve Švýcarsku  
DA Autoriseret repræsentant i Schweiz/importør i Schweiz  
DE Bevollmächtigter Vertreter in der Schweiz / Importeur in der Schweiz  
EL Εξουσιοδοτημένος εκπρόσωπος στην Ελβετία/ εισαγωγέας στην Ελβετία  
ES Representante autorizado en Suiza/Importador en Suiza  
FR Représentant autorisé en Suisse / Importateur en Suisse  
IT Rappresentante autorizzato in Svizzera/Importatore in Svizzera  
NL Geautoriseerde vertegenwoordiger in Zwitserland/importeur in Zwitserland  
NO Autorisert representant i Sveits / importer i Sveits  
PT Representante autorizado na Suíça/ Importador na Suíça  
SK Autorizovaný zástupca vo Švajčiarsku/dovozca vo Švajčiarsku  
SV Auktoriserad representant i Schweiz/importör i Schweiz  
TR İsviçre'deki Yetkili Temsilci/ İsviçre'deki İthalatçı  
ZH 瑞士授权代表/瑞士进口商

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	5
Anwendungsgebiete	6
Ενδείξεις χρήσης	7
Felhasználási javallatok	8
Indicazioni per l'uso	9
Bruksindikasjoner	10
Indicações de utilização	11
Indicaciones de uso	12
Indikationer för användning	13
Kullanım Endikasyonları	14
Indikace k použití	15
Návod na použitie	16
使用说明	17

**Indications for Use:**

The Posi-Stick® Introducer System is indicated for use in percutaneous introduction and placement of catheters and guidewires.

**Intended Patient Population and Group:**

For use on both adult and adolescent groups. The specific target patient population is dependent on the delivered device.

**Device Description:**

The Posi-Stick® Introducer consists of an inner stiffening cannula, dilator, and an outer dilator. The stiffening cannula accepts a .018" guidewire, and its design is to be inserted into the inner dilator. The device assembly may have a silicon or hydrophilic coating.

**Kit configurations may include the following devices:**

**Guidewire:** There are two configurations; Mandrel or Fully Coiled. Mandrel guidewire construction uses a tapered core wire which is soldered or welded to a coil at the distal end. Fully Coiled guidewire construction uses a tightly wound coil which surrounds a tapered core and a safety ribbon.

**Needle:** Introducer needles are composed of two components: a stainless steel cannula with an over molded hub. The introducer needle provides an access path into the target anatomy.

**Contraindications:**

Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel or target anatomy. There is increased risk of pneumothorax for the patient who has severe chronic lung disease.

**Potential Complications:**

The potential complications related to the use of the introducer include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis hematoma formation, hemothorax, vessel erosion, trauma to vessels, inadvertent arterial or venous puncture, local or systemic infection/sepsis.

**Precautions:**

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

**USA CAUTION**

**Federal (USA) law restricts this device to sale by or on the order of a physician.**

**Cautions:**

- This procedure should only be performed by physicians thoroughly trained in this procedure.
- Guidewires should be routinely inspected prior to use and discarded should any deformities be present in the guidewire.
- If resistance is met when advancing or withdrawing the guidewire or the introducer, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Because of the delicate and fragile nature of guidewires, extra care in handling must be taken.
- Do not expose to organic solvents, eg. alcohol. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

**Warnings:**

- Do not alter this device in any way.
- Do not reuse this device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- Do not attempt to straighten a wire that has been kinked or bent.
- Do not advance a guidewire that is kinked or becomes kinked or bent.
- Do not rotate the guidewire if significant resistance is felt.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.
- Do not soak in hydrophilic coated Posi-Stick® Introducer System in alcohol.
- To ensure lubricity, do not leave the hydrophilic coated Posi-Stick® Introducer System hydrated for long periods of time.
- Avoid abrasion of the coating of the hydrophilic coated Posi-Stick® Introducer System.
- During insertion avoid contact with bone, cartilage and scar tissue, which can damage the needle, guidewire or introducer

**USE STERILE TECHNIQUE.** A suggested procedure:

1. Peel open package and place contents on sterile field. Inspect introducer and accessories for defects. Do not use any defective devices.
2. To remove air, flush the introducer with normal saline.  
**For Hydrophilic Coated Posi-Stick® Introducer System only:** soak the introducer with saline for 15 to 20 sec. in order to completely wet and hydrate the surface. When hydrated (wetted with saline or blood), the introducer is very lubricious. Use sterile gauze soaked with saline to help with keeping the guidewire hydrated and to aid in its handling.
3. Prep skin and drape in area of anticipated puncture site as desired.
4. Insert 21ga trocar needle using sterile technique. The needle position should be verified.
5. The angle of the needle should be adjusted depending on the patient's build.
6. Remove trocar leaving needle cannula in place.
7. Advance the .018" guidewire through the introducer needle. Advance guidewire to required depth using ultrasound guidance, CT or fluoroscopy. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire location is suggested.
8. Hold guidewire in place and remove needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.  
**CAUTION:** Do not allow Guidewire to advance totally into patient.
9. Advance the Posi-Stick® Introducer over the .018" guidewire into the desired position.
10. Release and remove the stiffening cannula and the inner dilator, leaving the outer sheath and .018" guidewire in position.
11. Advance the .038" or .035" guidewire through the outer sheath (alongside the .018" guidewire) into the desired position.
12. Remove the .018" guidewire.
13. Remove the outer sheath leaving the .038" or .035" guidewire in place.
14. Utilize guidewire for placement of further diagnostic or interventional devices.

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**Emergo Europe**  
 Prinsessegracht 20  
 2514 AP The Hague  
 The Netherlands

CH REP

**MedEnvoy Switzerland**  
 Gotthardstrasse 28  
 6302 Zug  
 Switzerland

**UreSil, LLC**  
 5418 W. Touhy Ave  
 Skokie, IL 60077 USA  
 ☎ 1-800-538-7374 (U.S. only)  
 ☎ 1-847-982-0106

**Rx ONLY**

Caution: Federal law restricts this device to sale by or on the order of a (licensed healthcare practitioner)

STERILE EO

Ethylene Oxide Sterilized

Keep Dry

REF

Part Number

LOT

Lot Number

Manufacturer

NiTi

Nitinol

EC REP

Authorized Representative

NiTi [Au]

Nitinol Gold

NiTi [Pt]

Nitinol Platinum

Do Not Use If Package Is Damaged

Instructions for Use

Single Use Only

Caution

Not Made With Natural Rubber Latex

Non-Pyrogenic

Do Not Re-Sterilize

Store In A Dark Place

Use By Date

Manufacturing Date

Quantity

MD

Medical Device

MR Unsafe

Single Sterile Barrier system

Single Sterile Barrier System with Protective Packaging Outside

---

**DILATOR**

---

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	5
Anwendungsgebiete	6
Ενδείξεις χρήσης	8
Felhasználási javallatok	10
Indicazioni per l'uso	11
Bruksindikasjoner	12
Indicações de utilização	13
Indicaciones de uso	14
Indikationer för användning	15
Kullanım Endikasyonları	16
Indikace k použití	17
Návod na použitie	18
使用说明	19

---

## DILATOR

---

**Indications for Use:**

These dilators are used for the percutaneous introduction of guidewires into the peripheral vasculature.

**Device Description:**

The dilator introducer consists of a tipped radiopaque tube with a molded hub that accepts up to a .038 inch. guidewires

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**Contraindications:**

Use of the dilator is contraindicated if the patient has a known or suspected obstruction in the vessel. There is increased risk of pneumothorax for the patient who has severe chronic lung disease. Poor healing may result in the patient who has had irradiation to the anterior chest.

**Potential Complications:**

The potential complications related to the use of the dilator include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis, hematoma formation, hemothorax, vessel erosion, trauma to vessels, sepsis..

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**Cautions:**

- This procedure should only be performed by professionals thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire through the needle, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.
- Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

**Warnings:**

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

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1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture site as desired.
3. Distend the vessel of choice following standard hospital practice for venipuncture. The vein will be much easier to locate if the patient is well hydrated.
4. Insert needle into vessel. The needle position should be verified by observing venous blood return.
5. The angle of the needle should be adjusted depending on the patient's build: shallow in a thin person, deeper in a heavyset person.
6. Aspirate the puncture needle using the syringe.
7. Remove the syringe and insert soft tip of the guidewire through the introducer needle into the vessel. Advance guidewire to required depth. Leave an appropriate amount of guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire position is suggested.
8. Hold guidewire in place and remove introducer needle. Do not withdraw the guidewire back into the cannula as this may result in separation of the guidewire. The cannula should be removed first.  
**Caution: If the guidewire must be withdrawn while the needle is inserted, remove both the needle and wire as a unit to prevent the needle from damaging or shearing the guidewire.**
9. Thread the dilator over the guidewire.
10. Advance the dilator with a twisting motion over the guidewire and into the vessel, to dilate tissue in vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
11. Slowly withdraw and remove the dilator, while holding the guidewire in position and insert catheter or over guidewire.
12. Remove guidewire leaving catheter in position.

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