MD

UreSil®


EN

UreSil® TRU-CLOSE® Suction Drainage System

INTENDED USE: For collecting drainage fluids. The Suction Drainage Bags are intended to be used to collect fluids drained from abscesses, the hepatic system, and the kidneys and can also be used for post-operative wound drainage. This product is intended for short-term use (less than 30 days). Replace bag every 5 to 7 days or sooner if clinically indicated, for example, if malodorous or damaged.

INTENDED PATIENT POPULATION: This product is intended to be used with patients with placed drainage catheters or wound drains.

INTENDED USERS: This product is intended only for use as directed by a trained clinician. A nurse, caregiver, or patient may drain the fluids using the provided empty port (if present) and may compress the suction bellows to continue suction drainage as directed by a trained clinician.

NOTE: Graduation marks on bag not intended to measure – for reference only

CLINICAL BENEFITS: The Suction Drainage Bags are available in a variety of sizes and configurations. They can attach to any standard luer locking drainage catheter and can remove fluids using active suction or gravity. Additional clinical benefits:

- Simple and intuitive to operate.
- Compact and self-contained system allows full patient ambulation.
- During activation, body fluids are not aerosolized because all fluids are contained in the system.
- Closed system design prevents the chances of cross-contamination.
- Anti-reflux valve prevents liquids and air from backing up into the catheter.
- If the bellows fills, flow to the bag continues via gravity drainage as long as the bag is below the level of the drainage site.

PERFORMANCE CHARACTERISTIC: The Suction Drainage Bags are comprised of an inlet tube with a luer lock fitting for connection to a standard luer lock drainage catheter, internal anti-reflux valve, hydrophobic filter vents, suction bellows, and a Velcro belt system. Some versions of the bags have empty ports. Activation of the bellows initiates the suction drainage feature. Repeated activation of the bellows re-initiates the suction cycle. If the bellows fills, flow to the bag continues via gravity drainage as long as the bag is below the level of the drainage site. The bag is strong enough to withstand a buildup of air pressure within the bag. The bag is equipped with hydrophobic filter vents, which will vent collected air. Air can be manually forced out of the bag by gently squeezing the bag while the system is positioned vertically (housing above bag).

ACCESSORIES TO BE USED WITH THE DEVICE: The PDC-36 Percutaneous Drainage Catheter Connector Tubing is comprised of an inlet tube with a luer connector and is used for connecting a percutaneous catheter to a drainage bag or system. The TCY2-“Y” Connector Tubing Set is comprised of an inlet tube with a luer connector and a “Y” connector and is used for bilateral drainage with a gravity-drainage or suction-drainage system. UFA-LL is comprised of a tapered funnel and a female luer connector and is used to connect UreSil Ultra-flo Large Bore non-locking catheters (REF LBNL-series) to a standard luer.

INSTRUCTIONS FOR USE:

1. Close the empty port prior to use (if present).
2. Connect the drainage catheter to the luer lock as shown in Figure 2.
3. If altering the bag tube length, remove the tube from the connector at the top of the “suction bellows” housing, cut the tubing and re-attach.
4. Suspend the bag at a location below the drainage site.
5. Compress the bellows to aspirate fluid from the catheter as shown in Figure 3.
6. When bellows is full, repeat step 5 to reactivate suction.
7. Discard the entire system when the bag is full or when otherwise indicated.

CAUTION:

1. Use care if emptying bag to avoid contact with drained fluids.
2. Consult a healthcare professional if no fluid is draining into the bag or if fluid is leaking from the bag.

Bags without empty ports only: Do not attempt to empty the bag during use. The bag can be emptied prior to disposal if desired by cutting off the bottom corner of the bag.

CONTRAINDICATIONS: Currently no specific contraindications are known for the Suction Drainage Bags.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction. Sterile if package is unopened and undamaged. Do not use if the sterile package is damaged or is unintentionally opened before use.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil and TRU-CLOSE are registered trademarks of UreSil, LLC.

Figure 2


Figure 3

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