MD

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DA

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MD UreSil® "Y" Connector Tubing Set

EN

INTENDED USE: For connecting a percutaneous drainage catheter to a drainage bag or system. For bilateral drainage with a gravity drainage or suction drainage system. This product is intended for short-term use (less than 30 days).

INTENDED PATIENT POPULATION: This product is intended to be used with patients with placed drainage catheters.

INTENDED USERS: This product is intended only for use by trained clinicians.

"Y" Connector Tubing Set

CLINICAL BENEFITS: Work as an accessory in conjunction with a percutaneous drainage catheter and a Gravity Drainage bag or system to remove fluids using gravity or in conjunction with a percutaneous drainage catheter and a Suction Drainage bag or system to remove fluids using active suction or gravity.

PERFORMANCE CHARACTERISTIC: The TCY2 "Y" Connector Tubing Set is comprised of an inlet tube with a luer connector and a "Y" connector and is used to allow bilateral percutaneous catheter drainage into a single Gravity Drainage or Suction Drainage bag or system.

INTENDED FOR USE WITH: This accessory for a medical device is intended for use with bilateral percutaneous drainage catheters with standard luer lock and a Gravity Drainage or Suction Drainage bag or system.

INSTRUCTIONS FOR USE:

1. 2.

CONTRAINDICATIONS: Currently no specific contraindications are known for this accessory for a medical device.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction. Sterile if package is unopened and undamaged. Do not use if the sterile package is damaged or is unintentionally opened before use.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil is a registered trademark of UreSil, LLC.

REF EN Catalog number
CS Katalogové číslo
DA Katalognummer
DE Katalognummer
EL Αριθμός καταλόγου
ES Número de catálogo
FR Numéro de catalogue
IT Numero di catalogo
NL Catalogusnummer
PT Referência
SK Katalogové číslo
SV Katalognummer
TR Katalog numarası
ZH 产品目录号

LOT EN Batch code
CS Kód šarže
DA Batchkode
DE Chargencode
EL Κωδικός παρτίδας
ES Código de lote
FR Code du lot
IT Codice lotto
NL Coticdo do lote
PT Código do lote
SK Kód šarže
SV Batchkod
TR Parti kodu
ZH 批号

! EN Caution
CS Upozornění
DA Forsigtig
DE Achtung
EL Προσοχή
ES Precaución
FR Avertissement
IT Attenzione
NL Opgelet
PT Atenção
SK Upozornenie
SV Varning
TR Dikkat
ZH 注意

EN Country of manufacture/ Date of manufacture
CS Země výroby/ Datum výroby
DA Fremstillingsland/ Fremstillingsdato
DE Herstellungsland/ Herstellungsdatum
EL Χώρα κατασκευής/ Ημερομηνία κατασκευής
ES País de fabricación/ Fecha de fabricación
FR Pays de fabrication/ Date de fabrication
IT Paese di produzione/ Data di produzione
NL Land van productie/ Productiedatum
PT País de fabrico/ Data de fabrico
SK Krajina výroby/ Dátum výroby
TR Ambalaj hasarıyla ürünü kullanmayın ve kullanım talimatlarına başvurun/ Ürün kullanım tarihini ve kullanım talimatlarına başvurun/ Tilverkningsdato
SV Tillverkningsland/ Öretim tarihi
ZH 产地/生产日期

EN Use-by date
CS Datum použitelnosti
DA Sidste anvendelsesdato
DE Verfallsdatum
EL Ημερομηνία λήξης
ES Fecha de caducidad
FR Date de péremption
IT Data di scadenza
NL Houdbaarheidsdatum
PT Prazo de validade
SK Dátum spotreby
SV Sista förbrukningsdag
TR Son kullanma tarihi
ZH 有效期

EN Quantity
CS Množství
DA Antal
DE Menge
EL Ποσότητα
ES Cantidad
FR Quantité
IT Quantità
NL Hoeveelheid
PT Quantidade
SK Množstvo
SV Kvantitet
TR Miktar
ZH 数量


EN Not made with natural rubber latex
CS Není vyrobeno z latexu z přírodního kaučuku
DA Ikke fremstillet med naturgummlatex
DE Nicht aus Naturkautschuklatex hergestellt
EL Δεν κατασκευάζεται με λατέξ από φυσικό καουτσούκ
ES No fabricado con látex de caucho natural
FR Non fabriqué avec du latex de caoutchouc naturel
IT Non realizzato con lattice di gomma naturale
NL Niet gemaakt van natuurlijk rubberlatex
PT Não fabricado com látex de borracha natural
SK Neobsahuje prírodný latex
SV Inte tillverkad med naturlig gummlatex
TR Doğal kauçuk lateksten yapılmamıştır
ZH 并非天然橡胶材质

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EN Caution: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.
DE Achtung: Das Bundesgesetz der U.S.A. beschränkt den Verkauf und die Bestellung dieses Produkts auf Ärzte.