UreSii Aspirations-/skylladapter

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
INTENDED USE: For Catheter Aspiration and Flushing. This product is intended for short-term use (less than 30 days).


INTENDED PATIENT POPULATION: This product is intended to be used with patients with placed drainage catheters.

INTENDED USERS: This product is intended only for use by trained clinicians.

AFA-LL and DFA-LL for luer lock catheters:

AFA-LB for large bore catheters:

INSTRUCTIONS FOR USE:

- For catheter maintenance (aspiration and flushing), attach the appropriate adaptor to the catheter. The male threaded connector attaches to the catheter.
- Connect the other end of the adaptor to a collection bag. For the AFA-LB, remove the luer connector from the bag tubing and attach the bag tubing to the barb.
- Swab the silicone seal of the needleless port in accordance with facility protocol. Attach a male slip luer or male luer lock syringe to the needleless port by grasping the valve and firmly pushing and twisting the syringe onto the valve until tight.
- CAUTION:** Do not use needle to access the needleless port. For all adaptors other than the DFA-LL, flush or aspirate the catheter by clamping the tube with the supplied clamp and then aspirating/flushing as required. Unclamp the tube when complete and remove the syringe as required by twisting counter clock-wise until loose. Do not cap the valve. The device is closed when not engaged.
- The DFA-LL is provided with two clamps. When aspirating and flushing with the DFA-LL, close the back clamp and leave the front clamp open. For transferring fluid from the syringe to the bag, open the back clamp and close the front clamp. When aspiration and/or flushing is completed, remove the syringe and open both clamps. It is unnecessary to cap the needleless port.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil is a registered trademark of UreSil, LLC.

EN Not made with natural rubber latex
 CS Není vyrobeno z latexu z přírodního kaučuku
 DA Ikke fremstillet med naturgummlætex
 DE Nicht aus Naturkautschuklatex hergestellt
 EL Δεν κατασκευάζεται με Ατράξι από φυσικό καουτσόκ
 ES No fabricado con látex de caucho natural
 FR Non fabriqué avec du latex de caoutchouc naturel
 IT Non realizzato con lattice di gomma naturale
 NL Niet gemaakt van natuurlijk rubberlatex
 PT Não fabricado com látex de borracha natural
 SK Neobsahuje prírodný latex
 SV Inte tillverkad med naturlig gummlatex
 TR Doğal kauçuk latekten yapılmamıştır
 ZH 并非天然橡胶材质

EN Non-pyrogenic
 CS Apyrogenní
 DA Ikke-pyrogen
 DE Nicht pyrogen
 EL Μη πυρογόνο
 ES No contiene pirógenos
 FR Non pyrogène
 IT Non pirogenico
 NL Niet-pyrogeen
 PT Não pirogénico
 SK Nepyrogénne
 SV Iscke-pyrogen
 TR Pirojenik değildir
 ZH 无致热原

REF EN Catalog number
 CS Katalogové číslo
 DA Katalognummer
 DE Katalognummer
 EL Αριθμός καταλόγου
 ES Número de catálogo
 FR Numéro de catalogue
 IT Numero di catalogo
 NL Catalogusnummer
 PT Referência
 SK Katalogové číslo
 SV Katalognummer
 TR Katalog numarası
 ZH 产品目录号

LOT EN Batch code
 CS Kód šarže
 DA Batchkode
 DE Chargencode
 EL Κωδικός παρτίδας
 ES Código de lote
 FR Code du lot
 IT Codice lotto
 NL Batchcode
 PT Código do lote
 SK Kód šarže
 SV Batchkod
 TR Parti kodu
 ZH 批号

! EN Caution
 CS Upozornění
 DA Forsigtig
 DE Achtung
 EL Προσοχή
 ES Precaución
 FR Avertissement
 IT Attenzione
 NL Opgelet
 PT Cuidado
 SK Upozornenie
 SV Varning
 TR Dikkat
 ZH 注意

USA EN Country of manufacture/ Date of manufacture
 CS Země výroby/ Datum výroby
 DA Fremstillingsland/ Fremstillingsdato
 DE Herstellungsland/ Herstellungsdatum
 EL Χώρα κατασκευής/ Ημερομηνία κατασκευής
 ES País de fabricación/ Fecha de fabricación
 FR Pays de fabrication/ Date de fabrication
 IT Paese di produzione/ Data di produzione
 NL Land van productie/ Productiedatum
 PT País de fabric/ Data de fabrico
 SK Krajina výroby/ Dátum výroby
 TR Ambalaj hasarlıysa ürünü kullanmayın ve kullanım talimatlarına başvurun/ Ürünle ilgili kullanmayı ve kullanım talimatlarına başvurun/ İthalatçı
 SV Tillverkningsland/ Öretim tarihi
 ZH 生产地/ 生产日期

EN Use-by date
 CS Datum použitelnosti
 DA Sidste anvendelsesdato
 DE Verfallsdatum
 EL Ημερομηνία λήξης
 ES Fecha de caducidad
 FR Date de péremption
 IT Data di scadenza
 NL Houdbaarheidsdatum
 PT Prazo de validade
 SK Dátum spotreby
 SV Sista förbrukningsdag
 TR Son kullanma tarihi
 ZH 有效期

EN EN Quantity
 CS Množství
 DA Antal
 DE Menge
 EL Ποσότητα
 ES Cantidad
 FR Quantité
 IT Quantità
 NL Hoeveelheid
 PT Quantidade
 SK Množstvo
 SV Kvantitet
 TR Miktar
 ZH 数量

LATEX EN Not made with natural rubber latex
 CS Není vyrobeno z latexu z přírodního kaučuku
 DA Ikke fremstillet med naturgummlætex
 DE Nicht aus Naturkautschuklatex hergestellt
 EL Δεν κατασκευάζεται με Ατράξι από φυσικό καουτσόκ
 ES No fabricado con látex de caucho natural
 FR Non fabriqué avec du latex de caoutchouc naturel
 IT Non realizzato con lattice di gomma naturale
 NL Niet gemaakt van natuurlijk rubberlatex
 PT Não fabricado com látex de borracha natural
 SK Neobsahuje prírodný latex
 SV Inte tillverkad med naturlig gummlatex
 TR Doğal kauçuk latekten yapılmamıştır
 ZH 并非天然橡胶材质

XX EN Non-pyrogenic
 CS Apyrogenní
 DA Ikke-pyrogen
 DE Nicht pyrogen
 EL Μη πυρογόνο
 ES No contiene pirógenos
 FR Non pyrogène
 IT Non pirogenico
 NL Niet-pyrogeen
 PT Não pirogénico
 SK Nepyrogénne
 SV Iscke-pyrogen
 TR Pirojenik değildir
 ZH 无致热原

EC REP EMERGO EUROPE, Westervoortseijk 60, 6827 AT Arnhem, The Netherlands

EN Authorized representative in the European Community / European Union
 CS Způsobilý zástupce v Evropské společnosti / Evropské unii
 DA Autoriseret repræsentant i Det Europæiske Fællesskab/Den Europæiske Union
 DE Bevollmächtigter in der Europäischen Gemeinschaft / Europäischen Union
 EL Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα / Ευρωπαϊκή Ένωση
 ES Representante autorizado en la Comunidad Europea / Unión Europea
 FR Représentant autorisé dans la Communauté européenne et l'Union européenne
 IT Rappresentante autorizzato nella Comunità europea/Unione Europea
 NL Geautoriseerde vertegenwoordiger in de Europese Gemeenschap / Europese Unie
 PT Representante autorizado na Comunidade Europeia / União Europeia
 SK Spôsobný zástupca v Európskom spoločenstve alebo Európskej únii
 SV Auktoriserad representant inom EU/EES
 TR Avrupa Topluluğu / Avrupa Birliği/deki yetkili temsilci
 ZH 欧洲共同体/欧盟的授权代表

EC MedEnvoy, Prinses Margrietplantsoen 33 Suite 123, 2595 AM, The Hague, The Netherlands

EN Importer in the European Community / European Union
 CS Dovezoc v Evropské společnosti / Evropské unii
 DA Importør i Det Europæiske Fællesskab/Den Europæiske Union
 DE Importeur in der Europäischen Gemeinschaft / Europäischen Union
 EL Εισαγωγής στην Ευρωπαϊκή Κοινότητα / Ευρωπαϊκή Ένωση
 ES Importador en la Comunidad Europea / Unión Europea
 FR Importateur dans la Communauté européenne et l'Union européenne
 IT Importatore nella Comunità europea/Unione Europea
 NL Importeur in de Europese Gemeenschap / Europese Unie
 PT Importador na Comunidade Europeia / União Europeia
 SK Dovezoc v Európskom spoločenstve alebo Európskej únii
 SV Importör inom EU/EES
 TR Avrupa Topluluğu / Avrupa Birliği/ndeki İthalatçı
 ZH 进口商欧洲共同体/欧盟的

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