
UreSii TRU-FIX Kateterfiksationssæt

DA

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
INTENDED USE: For externally securing percutaneous catheters. This product is intended for short-term use (less than 30 days).

INTENDED USERS: This product is intended only for use by trained clinicians.

INSTRUCTIONS FOR USE:

1. Clean and dry the catheter exit site.
2. Position the disk onto the catheter via the slit.
3. Slide the disk onto the skin.
4. Close the clamp completely.
5. Apply the adhesive dressing after removing the release paper.
6. The adhesive patch is expected to adhere for at least 7 days. If replacement of patch is warranted or desired, remove carefully by peeling away from skin. Clean and dry the catheter exit site and place a new patch (REF TFD-R).

CONTRAINDICATIONS:


1. Not intended for vascular use.
2. Do not use on patients with known tape or adhesive allergies.

WARNING: The reuse of this single-use device can lead to patient infection and/or device malfunction.

SAFE DISPOSAL: Dispose of used device in container marked for biohazard (i.e. contaminated with potentially infectious substances of human origin).

REPORTING SERIOUS INCIDENTS: Any serious incident that has occurred in relation to the device should be reported to UreSil as the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

UreSil is a registered trademark of UreSil, LLC. TRU-FIX is a trademark of UreSil, LLC. U.S. Patent no. 5,833,666.

EN Not made with natural rubber latex
CS Není vyrobeno z latexu z přírodního kaučuku
DA Ikke fremstillet med naturgummilætex
DE Nicht aus Naturkautschuklætex hergestellt
EL Δεν κατασκευάζεται με latex από φυσικό καουτσόκ
ES No fabricado con látex de caucho natural
FI Ei ole valmistettu luonnonkumilætekistä
FR Non fabriqué avec du latex de caoutchouc naturel
IT Non realizzato con lattice di gomma naturale
NL Niet gemaakt van natuurlijk rubberlatex
NO Skal ikke brukes med naturlig gummlætex
PT Não fabricado com látex de borracha natural
SK Neobsahuje prírodný latex
SV Inte tillverkad med naturlig gummlætex
TR Doğal kauçuk lâtekten yapılmamıştır
ZH 并非天然乳胶材质

REF EN Catalog number
CS Katalogové číslo
DA Katalognummer
DE Katalognummer
EL Αριθμός καταλόγου
ES Número de catálogo
FI Luettelonumero
FR Numéro de catalogue
IT Numero di catalogo
NL Catalogusnummer
NO Katalognummer
PT Referência
SK Katalógové číslo
SV Katalognummer
TR Katalog numarası
ZH 产品目录号

LOT EN Batch code
CS Kód šarže
DA Batchkode
DE Chargencode
EL Κωδικός παρτίδας
ES Código de lote
FI Eräkkoodi
FR Code du lot
IT Codice lotto
NL Batchcode
NO Partikode
PT Código do lote
SK Kód šarže
SV Batchkod
TR Parti kodu
ZH 批号

! EN Caution
CS Upozornění
DA Forsigtig
DE Achtung
EL Προσοχή
ES Precaución
FI Huomio
FR Avertissement
IT Attenzione
NL Opgelet
NO Advarsel
PT Atenção
SK Upozornenie
SV Varning
TR Dikkat
ZH 注意

USA EN Country of manufacture/ Date of manufacture
CS Země výroby/ Datum výroby
DA Fremstillingsland/ Fremstillingsdato
DE Herstellungsland/ Herstellungsdatum
EL Χώρα κατασκευής/ Ημερομηνία κατασκευής
ES País de fabricación/ Fecha de fabricación
FI Valmistusmaa/ Valmistuspäivä
FR Pays de fabrication/ Date de fabrication
IT Paese di produzione/ Data di produzione
NL Land van productie/ Productiedatum
NO Produksjonsland/ Produksjonsdato
PT País de fabrico/ Data de fabrico
SK Krajina výroby/ Datum výroby
TR Ambalaj hasarlıysa ürünü kullanmayın ve kullanım talimatlarınıza başvurun/ Tillverkningsdag
SV Tillverkningsland/ Öretim tarihi
ZH 生产地/ 生产日期

EN Use-by date
CS Datum použitelnosti
DA Sidste anvendelsesdato
DE Verfallsdatum
EL Ημερομηνία λήξης
ES Fecha de caducidad
FI Viimeinen käyttöpäivä
FR Date de péremption
IT Data di scadenza
NL Houdbaarheidsdatum
NO Bruk før-dato
PT Prazo de validade
SK Dátum spotreby
SV Sista förbrukningsdag
TR Son kullanma tarihi
ZH 有效期

EN Quantity
CS Množství
DA Antal
DE Menge
EL Ποσότητα
ES Cantidad
FI Määrä
FR Quantité
IT Quantità
NL Hoeveelheid
NO Antall
PT Quantidade
SK Množstvo
SV Miktat
TR Miktar
ZH 数量

LATEX EN Not made with natural rubber latex
CS Není vyrobeno z latexu z přírodního kaučuku
DA Ikke fremstillet med naturgummilætex
DE Nicht aus Naturkautschuklætex hergestellt
EL Δεν κατασκευάζεται με latex από φυσικό καουτσόκ
ES No fabricado con látex de caucho natural
FI Ei ole valmistettu luonnonkumilætekistä
FR Non fabriqué avec du latex de caoutchouc naturel
IT Non realizzato con lattice di gomma naturale
NL Niet gemaakt van natuurlijk rubberlatex
NO Skal ikke brukes med naturlig gummlætex
PT Não fabricado com látex de borracha natural
SK Neobsahuje prírodný latex
SV Inte tillverkad med naturlig gummlætex
TR Doğal kauçuk lâtekten yapılmamıştır
ZH 并非天然乳胶材质

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