




Global Product Catalog

CT, Interventional Radiology, Ultrasound
and Emergency Department





Origin[®]

General Purpose GP[™] Drainage Catheter with locking pigtail



- Origin Locking Catheter features a clear hub for visualization, tapered tip, large holes, and depth marks that can indicate catheter position.
- Every Origin General Purpose Catheter (GP) includes the following: Catheter, Metal Cannula, Flexible Cannula, Trocar and Hydrophilic Coating. The 6F catheter also includes an obturator.
- Origin Catheters are available in 6F through 16F diameters.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Size/Length	Units/Box
GPL2-0625H	6 F / 18 cm	5
GPL2-0830H	8 F / 21 cm	5
GPL2-1030H	10 F / 21 cm	5
GPL2-1230H	12 F / 20 cm	5
GPL2-1430H	14 F / 21 cm	5
GPL2-1630H	16 F / 21 cm	5

Origin[®] Mini-Pigtail[®]

Drainage Catheter with locking pigtail



- UreSil's "MiniPig" (MP) Catheter is very similar to the General Purpose Catheter but is designed with a reduced size pigtail that presents a smaller footprint within an abscess, and is generally used for abscesses near critical structures.
- Every MiniPig General Purpose Catheter includes the following: Catheter, Metal Cannula, Flexible Cannula, Trocar and Hydrophilic Coating. The 6F catheter also includes an obturator.
- MiniPig Catheters are available in 6F, 8F and 10F diameters.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Size/Length	Units/Box
MPL2-0625H	6 F / 20 cm	5
MPL2-0830H	8 F / 25 cm	5
MPL2-1030H	10 F / 25 cm	5

Origin® Short Mini-Pigtail®

Drainage Catheter with locking pigtail



- UreSil's Short MiniPig Catheter (KDL) has a working length of only 12cm.
- This catheter is suitable for the drainage of abscesses closer to the surface or in patients with smaller anatomy.
- Every KDL Catheter includes the following: Catheter, Metal Cannula, Trocar and Hydrophilic Coating
- Short MiniPig Catheters are available in 6F, 8F, and 10F diameters
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Size/Length	Units/Box
KDL2-0617H	6 F / 12 cm	5
KDL2-0817H	8 F / 12 cm	5
KDL2-1017H	10 F / 12 cm	5

Origin[®] General-Purpose

Drainage Catheter with “J” tip



- UreSil's Non-Locking General Purpose Catheter (GPJ) features a tapered tip, large drain holes, and depth marks that can indicate catheter position.
- Every GPJ Catheter includes the following: Catheter, Metal Cannula, Trocar and Hydrophilic Coating.
- GPJ Catheters are available in 6F through 14F diameters.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Size/Length	Units/Box
GPJ-0620H	6 F / 15 cm	5
GPJ-0825H	8 F / 19 cm	5
GPJ-1025H	10 F / 19 cm	5
GPJ-1225H	12 F / 19 cm	5
GPJ-1425H	14 F / 20 cm	5



Origin[®] NEPHRO-CATH[™]

Neophrostomy Drainage Catheter
with locking pigtail



- UreSil's Nephrostomy Catheter (NCL) features a clear hub for visualization, tapered tip, large drain holes, and depth marks that can indicate catheter position. The "35" version of this catheter is approximately 5 cm longer than General Purpose Catheters to allow for placement within the kidney for drainage.
- UreSil has a longer "50" version of this catheter to accommodate larger patients where increased length is necessary.
- NCL Catheter is available in 8F-14F diameters and includes the following: Catheter, Metal Cannula, Flexible Cannula and Septum cap (no trocar).
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Size/Length	Units/Box
NCL2-0835H	8 F / 28 cm	5
NCL2-1035H	10 F / 28 cm	5
NCL2-1235H	12 F / 27 cm	5
NCL2-1435H	14 F / 28 cm	5

Product Code	Size/Length	Units/Box
NCL2-0850HB	8 F / 41 cm	5
NCL2-1050HB	10 F / 41 cm	5
NCL2-1250HB	12 F / 40 cm	5
NCL2-1450HB	14 F / 40 cm	5



Origin[®] CHOLE-CATH[®]

Biliary Drainage Catheter
with locking pigtail



- UreSil offers a wide array of Biliary Catheters (CCL) – 8F to 16F diameters.
- All Biliary Catheters come with a radiopaque (RO) marker located just proximal to the proximal hole to facilitate proper placement under imaging guidance.
- All regular (HB) Biliary Catheters are available in 8F to 16F diameters, and have holes up the shaft approximately 9.5 cm from the curled pigtail.
- All extended hole Biliary Catheters (HBX) are available in 8F-12F diameters with additional holes extending up the shaft approximately 16 cm from the curled pigtail.
- Every Biliary Catheter includes the following: Catheter, Metal Cannula, Flexible Cannula, and Septum Cap
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Size/Length	Units/Box
CCL2-0850HB	8 F / 41 cm	5
CCL2-1050HB	10 F / 41 cm	5
CCL2-1250HB	12 F / 40 cm	5
CCL2-1450HB	14 F / 40 cm	5
CCL2-1650HB	16 F / 40 cm	5

Product Code	Size/Length	Units/Box
CCL2-0850HBX	8 F / 41 cm	5
CCL2-1050HBX	10 F / 41 cm	5
CCL2-1250HBX	12 F / 40 cm	5



Origin[®]CHOLE-CATH[®]

Biliary Drainage Catheter
with locking pigtail



Origin CHOLE-CATH Biliary Drainage Catheter Kit with RO Marker

UreSil's Biliary Kits include items needed to perform a biliary procedure.

Product Code	Size/Length	Units/Box
CCL2-0850HB-K	8 F / 41 cm	5
CCL2-1050HB-K	10 F / 41 cm	5

Origin NEPHRO-CATH[™]

Nephrostomy Drainage Catheter with
locking pigtail



Origin NEPHRO-CATH Nephrostomy Drainage Catheter Kit

UreSil's Nephrostomy Kits include items needed to perform a nephrostomy procedure.

Product Code	Size/Length	Units/Box
NCL2-0835HB-K	8 F / 28 cm	5
NCL2-1035HB-K	10 F / 28 cm	5

The kits are available in 8F and 10F diameters and include the following:

UreSil Biliary Catheter (CCL2-0850HB or CCL2-1050HB) or UreSil Nephrostomy Catheter (NCL2-0835H or NCL2-1035H) comprised of:

- Catheter
- Metal Cannula
- Flexible Cannula
- Septum Cap


Posi-Stick[®] Introducer Kit comprised of:

- Introducer Set with Radiopaque Tip (Coaxial 4F/6F dilator set with cannula)
- 21 GA needle with Echogenic Tip
- .018", 70cm long Nitinol Guidewire with Floppy Radiopaque Tip
- .038", 150cm long Guidewire


Dilator Pack comprised of:

- 8F, 10F, and 12F Dilators

Procedure packs are comprised of CE-marked sterile devices and are available for sale in the USA and in countries where the product has received clearance for sale where applicable.



THORA-VENT®



Red Signal Diaphragm

- Fluctuations demonstrate entry of the catheter into the pleural cavity with pressures greater than 6mm Hg and less than - 6mm Hg
- Continues to fluctuate with respiration (as the intrapleural pressure varies from > +6mmHg to < -6mmHg) until the pneumothorax is resolved as confirmed by x-ray

- The Thora-Vent is designed for the treatment of spontaneous, traumatic or iatrogenic simple pneumothorax.
- Thora-Vent is unique in that it allows for active monitoring of the patient's pneumothorax during treatment and provides for maximum ambulation as well.
- Available in multiple French sizes and lengths, Thora-Vent is an excellent option for the resolution of simple pneumothorax.
- Every kit includes the following: Thora-Vent with Trocar, Aspiration Cannula, Suction Tubing Set, 50cc Syringe, Halstead Forceps, Vent Occlusion Plug, Gauze Sponges, Protected Disposable Scalpel, and Fenestrated Drape.
- *Thoracic Vent Procedure Tray, TVST, and TVOCP are CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.
- **Available over the wire insertion accessory (REF OTW11-13) is available for sale in the USA and in countries where the product has received clearance where applicable. Product is not CE marked and is therefore not available where CE mark is required.

*THORA-VENT Procedure Tray

Product Code	Size/Length	Units/Box
TV11-10	11 F / 10 cm	1
TV11-13	11 F / 13 cm	1
TV13-10	13 F / 10 cm	1
TV13-13	13 F / 13 cm	1


THORA-VENT Accessories

Product Code	Description	Units/Box
*TVST	Thora-Vent Suction Tubing	5
*TVOCP	Vent Occlusion Plug	5
**OTW11-13	Over-the wire insertion Cannula for TV11-13	5



Evolution®

Evacuated Suction Bottles



ESB-1000SP Non-Sterile

- The Evolution line of evacuated suction bottles is used for paracentesis, thoracentesis and other drainage procedures.
- The Evolution Bottle is not made from glass.
- The available male-to-male luer connecting tubing allows for flow without a spike.
- Evolution is not made with natural rubber latex.
- These bottles are available in a non-sterile format with sterile tubing sold separately or packed sterile with tubing already connected.
- The unique roller clamp design provides an easy way to control the flow of fluids if desired
- Available with spikeable pour spout.
- Product is available for sale in the USA and in countries where the product has received clearance for sale where applicable. Product is not CE marked and is therefore not available where CE mark required for sale.



ESB-1000FS-T Sterile with Drain Line

Product Code	Size/Length	Units/Box
ESB-1000SP	1000 ml Evacuated Bottle w/ Luer Connection Only, Non-Sterile	12
ESB-1000FS-T	Sterile 1000 ml Evacuated Bottle w/ Drain Line (120 cm)	10
ESB-1001SP	1000 ml Evacuated Bottle w/ Luer Connection and spikeable, pour spout, Non-Sterile	12
ESB-ST	Sterile Tubing Only (120 cm)	10
ESB-ST1	Sterile Tubing with spike (120 cm)	10




Posi-Stick®

Introducer Sets



One Step™ Introducer Set



Product Code	Diagnostic Needle (21 GA)	*POSI-STICK® Introducer Set	Nilinol Guidewire (0.018" 70 cm length)	Guidewire (0.038" 150 cm length)	Guidewire (0.035" 150 cm length)	Dilator (5 Fr) with locking trocar	Guidewire (0.038" 100 cm length)	Units per Box
PS2-38	✓	✓	✓	✓				5
PS2-35	✓	✓	✓			✓		5
PS2-18*			✓					10
PS2-00	✓	✓	✓					10
PS2-15							✓	5

* Replacement 0.018" (.46 mm) Guidewire for PS2-35 & PS2-38

* Replacement 0.018\"/>

- PS2-38, PS2-35, PS2-18, and PS2-00 are CE-marked.



Manufactured by Galt Medical Corp, 2220 Merritt Drive,
Garland, TX 75041 USA **CE**₁₆₃₉



Authorized Representative in the European Community:
Emergo Europe, Westervoortsedijk 60, 6827 AT Arnhem,
The Netherlands

- PS2-1S: Product is available for sale in the USA and in countries where the product has received clearance for sale where applicable. Product is not CE marked and is therefore not available where CE mark required for sale.

"POSI-STICK" Set Configurations



TRU-CLOSE®

TCS Suction Drainage Bag

- UreSil's TRU-CLOSE TCS Suction Drainage Bag will attach to any luer locking drainage catheter and can remove fluids using active suction or gravity.
- Available in a variety of sizes and configurations, every TRU-CLOSE bag has an anti-reflux valve, hydrophobic filter vents, provides active suction and is simple to operate.
- TCS Bags are available in 300ml, 500ml and 1,000ml configurations.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Description	Units/Box
TCS500	500 ml w/o Drainage Port	10
TCS1000	1000 ml w/o Drainage Port	10
TCS500D	500 ml with Drainage Port	10
TCS300D	300 ml with Drainage Port	10
TCS500DS	Low profile bellows	10
	500 ml with Drainage Port	
	Low profile bellows	

TCS500DS



TCS1000



TCS300D



Optimus®

Gravity Drainage Bag




Product Code	Description	Units/Box
TC600	Gravity Drainage System with 600 ml drainable collection bag, elastic hook-and-loop belt and 20" stretchable tubing with luer lock connector	20
TC600L	Gravity Drainage System with 600 ml drainable collection bag, elastic hook-and-loop belt and 36" stretchable tubing with luer lock connector	20

- UreSil's Optimus Gravity Drainage Bags are made from plastic film and provide collection for interventional drainage.
- The TC600 comes with a 600ml drainable collection bag, 20" of stretchable, resizable tubing, a luer lock connector, a loop belt and an elastic hook.
- The TC600L is identical, but comes equipped with 36" of stretchable, resizable tubing with a luer lock connector.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.


Aspiration Flush Adapters

AFA-LL



- The Aspiration Flush Adapter with Luer Lock Connections replaces the flow limiting three way stop-cock with something easy to use.
- The AFA-LL can be used with any luer lock bag and any luer lock catheter.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

AFA-LB




- The AFA-LB is identical to the device provided in our LBL Catheter.
- The AFA-LB allows connection of the UreSil LBL catheter to any standard drainage bag.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Units/Box
AFA-LL	5


Product Code	Units/Box
AFA-LB	5

DFA-LL



- Similar to the AFA-LL but includes an additional clamp that can be used to allow flushing towards the catheter or the bag and the catheter can be clamped to ensure flushing towards the bag.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

UFA-LL




- This device provides a way to adapt funneled catheter connections to a luer lock drainage bag or system.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Units/Box
DFA-LL	5

Product Code	Units/Box
UFA-LL	5


PDC-36



- PDC-36 multipurpose tubing is 36" in length. Common uses include connecting a basic urinary drain bag or water seal device to any standard drainage catheter by utilizing a luer connection.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Units/Box
PDC-36	10

TCY2




- The TCY2 tubing is used to offer bilateral drainage to a single bag.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	Units/Box
TCY2	5

TRU-FIX™

Fixation Devices



- The TRU-FIX Catheter Fixation Device provides catheter fixation.
- Available in two sizes to suit your needs, this product includes: Fixation Disk with Securing Clamp and Adhesive Dressing.
- Replacement Adhesive Dressing also available.
- Product is CE-marked and available for sale in the USA and in countries where the product has received clearance for sale where applicable.

Product Code	French Size	Units/Box
TFK-08	8F	5
TFK-14	14F	5

Replacement Adhesive Dressing Product Code	Units/Box
TFD-R	20

Refer to specific product pages for product's regulatory status and availability for sale. Where product is referenced as CE-marked, the following details are applicable, except where specifically noted:



Manufactured by UreSil, LLC, 5418 W Touhy Ave, Skokie, IL 60077 USA

Authorized Representative in the European Community:



Emergo Europe, Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands



THORA-VENT, Mini-Pigtail, UreSil, CHOLE-CATH, Optimus, Evolution, TRU-CLOSE, and POSI-STICK are registered trademarks of UreSil, LLC. **GP, NEPHRO-CATH, TRU-FIX, and One Step** are trademarks of UreSil, LLC.

Products covered by one or more of the following U.S. patents: 8,177,773; 9,079,006B1

Doc. no. MK-UPC Rev 8/10/2026 © 2026 UreSil, LLC

Terms are Net 30 Days. All shipments are F.O.B. Skokie, Illinois.

Prices are subject to change without notice.

RETURNED GOODS:

Approval to return merchandise is required for all returns. Contact UreSil Customer Service for a Return Goods Authorization Number at 800-538-7374. The request for returns should include the reason, catalog number, lot number and UreSil Invoice number. All items must be returned unopened, in full boxes only and are subject to a 20% restocking charge. Returns will only be accepted within ninety days of date of invoice. Material that has been opened, damaged or tampered with is not returnable for credit.





